

ClearKlens IPA

Sterile isopropyl alcohol for disinfection in cleanrooms

Technical data

Appearance	Clear liquid free from insolubles
pH	6-9
Density (20°C)	0.865-0.885
Alcohol proportion	70% +/- 1% v/v

The above data is typical of normal production and should not be taken as a specification.

Safe handling and storage information

Store in original closed containers, away from sunlight and extremes of temperature. Full guidance on the handling and disposal of this product is provided in a separate Material Safety Data Sheet.

Product compatibility

Applications solutions, when used as directed, will not affect materials normally encountered in the Pharma and Cosmetics industries.

Test method

Each lot is provided with:

Certificate of Analysis, available also on www.clearklens.com

Certificate of Sterility

Certificate of Irradiation