

# Autoclave Test Report

## PRODUCT

Pre-box

## Composition

100% Polypropylene



## TEST CONDITIONS

### Autoclave

Varioklav 135 S, steam sterilizer, single pre-vacuum, usable volume 135 litres

### Test Parameter

121°C for 20 minutes

Pre-vacuum: 20 kPa

Cooling-down of product between autoclaving cycles

## RESULT

After 80 autoclaving cycles the Pre-box does not show any damages and fulfills FHP specification requirements.



Photo: Pre-box after 80 autoclaving cycles

## ADDITIONAL COMMENTS

- This test was done in the laboratory of Freudenberg HP KG, Regional Technical Centre.
- The product was not used for cleaning between autoclaving cycles.
- Different conditions at the customer can lead to a different result.

**Date: 2011/07/18**

Distribuido por:



García Martín 19  
28224 Pozuelo, Madrid  
T: 917 091 540 | [ventas@anade.com](mailto:ventas@anade.com)  
[anade.com](http://anade.com)